


## Sphere Characterization Image Analysis Report 227

### Sample Description

Image of spheres from a SEM in BMP format.

### Purpose of Analysis


Demonstrate the ability of the Clemex Vision image analyzer to discriminate and measure the spheres and their pores.



**Figure 1:** a) Original image. b) The original image was modified by a gray transformation to outline the spheres (Sobel).

### Procedure


Some binary operations were used to complete spheres, remove artifacts and separate spheres (Closing, Chord Size, Opening, Separate).



**Figure 2:** a) The sphere outlines were binarized into red bitplane by Gray Thresholding. b) Final detection of spheres (red) and pores (blue) overlaid against the modified gray image.

### Results Summary

	Diameter ( $\mu\text{m}$ )	Sphericity
Minimum:	45.5	0.87
Maximum:	180.7	0.94
Mean:	111.3	0.91
Standard Deviation:	59.7	0.03
Count:	8	
Pores (area %)	3.84	



**Figure 3:** a) Circular Diameter distribution. b) Area of pores compared to the total area covered by the spheres (as a percentage).

### Equipment

#### Image Analysis

System:  
Source :

Clemex Vision SE  
Image from a SEM in BMP format.