

Image Analysis Report # 796

Bubble Size Measurement
for
Nova Chemicals

Samples submitted

2 images of bubbles in glycerin/mineral oil emulsion

General analysis request

Measure diameter of bubbles


Analysis procedure

Steps


1. The bubbles are detected based on their contrast with the background.
2. Touching bubbles are separated based on changes in convexity along their surface.
3. Small objects, considered as noise, are removed (3 x 3 pixels)
4. bubbles intersected by the field of view are excluded from the analysis.
5. Diameter and volume are measured for all bubbles.
6. Images loaded in groups are analyzed together.
7. Charts and statistics are generated automatically.


Charts and statistics

By count



by volume





by count

Minimum: 58.4 μm
 Maximum: 3333.1 μm
 Mean: 189.2 μm
 Std Dev.: 226 μm
 Sum: 780736 μm
 Count: 4126
 Under: 0
 Over: 0
 Accepted: 100 %
 Field Count: 2
 Field Area: 4.74E+08 μm^2
 Total Area: 9.48E+08 μm^2
 D10: 68 μm
 D50: 114.5 μm
 D90: 369.7 μm

by volume

Minimum: 58.4 μm
 Maximum: 3333.1 μm
 Mean: 1474.2 μm
 Std Dev.: 822.2 μm
 Sum: 2.70E+14 μm
 Count: 4126
 Under: 0
 Over: 0
 Accepted: 100 %
 Field Count: 2
 Field Area: 4.74E+08 μm^2
 Total Area: 9.48E+08 μm^2
 D10: 592.5 μm
 D50: 1333.6 μm
 D90: 2327.9 μm

Pictures of detections

